

Nursing Training Simulator

Equipment: Quest 3/Quest 2

Set up:

1. Turn on the headset and create a new boundary. The playspace should be standing and allow for full range of arm motion.
2. Open the App Library then locate and select Nursing Sim.

Controls/Navigation:

- Trigger to teleport
- Right trigger to select
- Grip to leave/start evaluation and grab

Gameplay:

1. The app will open with a “welcome” screen in front of the player, a wall behind the screen with a window peering into a nurse and patient. The screen will ask if the player wants to do a tutorial or to start the experience
 1. The tutorial doesn’t work and isn’t fully accurate, so just continue on to the experience
2. The player will then be asked to put in their User ID. After clicking enter, they’ll be taken to a screen that has a menu bar underneath the big window that says “Welcome to Nursing VR” and has an “About” and “Credits” button
 1. The menu bar has the time and the headset’s battery percentage on the left and the options “about”, “modules”, “results”, “tutorial”, and “settings”
3. To begin the experience, the player should select “Modules” to see the modules available (will be a mix of purple and blue colors if available). After selecting a module (listed on the left side of the screen), submodules should appear in the middle of the screen. Select an available submodule then select “begin simulation” on the right side of the screen to start the experience
4. The player will then follow the directions provided on the screen before going to the door and “leaving” to start their evaluation to see how much they remember from the practice.
 1. Many functions during practice and evaluation don’t work, so the player might not be able to get through the whole module.

Revision #1

Created 2026-10-07 18:25:52 UTC by Febronia Farowez

Updated 2026-10-07 18:26:38 UTC by Febronia Farowez